Peiying®

Alkomat AT551
MIE2124

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User's manual **EN**

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Manual de utilizare **RO**

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Read this instruction manual before use and keep it for future reference.

Introduction

- This device is intended for measuring alcohol concentration in human breath. Measurements obtained by this device are used in the diagnosis of alcohol intoxication.
- The device is for personal use, must not be used for medical or diagnostic purposes!

Preparation for Correct Measurement


- Before testing, wait at least 15 minutes after your last drinking.
- Only use the device in a temperature 10°C - 40°C (50°F - 104°F).


Product description

1. LCD screen
2. Mouthpiece
3. Battery cover
4. Sticker
5. Unit
6. Test button
7. Clip

Power supply


- To change battery, slide back the battery cover on the back of the device.
- Insert 1x AAA battery, making sure the polarization is correct.
- If the power is too low, the unit will display a battery symbol LP; the device may not work properly. Replace the battery.

Operation

- Insert one of the mouthpieces into the ,blow in' hole.
- Press and hold the test button until the device powers on.

- Wait for countdown to zero.
- Once the countdown is over, the screen will display flashing C letter. Take a deep breath, blow steady into the mouthpiece until hearing a CLICK sound (about blowing 6 seconds).

- The screen displays 0; wait for the test results.

- Your test result will be displayed on the screen for 20 seconds.

- The LCD screen will change the color if the test result is above alarm point.
- It will show error notification Er-b when incorrect blowing (refer to the chapter of the “Troubleshooting” for details). The device will be turned off after indicated error.

Troubleshooting

Error message/Problem	Description	Remarks
No display	<ul style="list-style-type: none"> ·Press and hold the TEST button until power on ·Check the battery; use alkaline battery only 	
Auto off	<ul style="list-style-type: none"> ·Low battery ·Over the period of test stage 	
Err-b	<ul style="list-style-type: none"> ·Blowing unsteady or interrupt ·Blowing less than 6 seconds ·Blowing flow too weak ·No blow during BLOW mode 	
LP	<ul style="list-style-type: none"> ·Low battery indication ·Need to change new battery ·Use alkaline battery only 	
HI	<ul style="list-style-type: none"> ·The blood alcohol content above the display range 	

Important notes

- Wait 15 minutes to test since your last drinking.
- The blowing time must not be below required time (around 6 seconds).
- The screen may display error if there's no breath blow, it has been interrupted or too short.
- If the device changes temperature suddenly, wait some minutes to use.
- Device should be used in temperature 10°C ~ 40°C (50°F ~ 104°F).
- Recalibration interval: 12 months or 200 tests
- Never store the device in a corrosive gas room.
- Do not expose the device to direct sun light nor sources of heat.
- Avoid dropping or subjecting the device to a strong shock.
- Selective unit of measure:

To change unit of measurement

- Press and hold TEST button.
- The device will begin countdown.
- Once the countdown is over, the screen will display flashing C letter (do not blow into the mouthpiece).
- Next, the screen will display error notification Er-b - press and hold TEST button.
- The screen will display UnIt. Units will be displayed on the right side of the screen. To change units repeatedly press TEST button (available units: %, mg/L, ‰).
- In order to save, leave the screen displaying desired unit and wait until device emits sound signal and turns off – change has been successfully saved; OR press and hold the TEST button again – the device will enter selective alarm point function.
- Selective alarm point:

To change alarm point

- press and hold TEST button.
- The device will begin countdown.
- Once the countdown is over, the screen will

- display flashing C letter (do not blow into the mouthpiece).
- Next, the screen will display error notification Er-b - press and hold TEST button.
 - The screen will display Unlt. Press and hold the TEST button again. The screen will display ALAr.
 - To change alarm point, repeatedly press TEST button; the device will display alarm point value (available alarm points for chosen units: ‰ - 0,020, 0,030, 0,050, 0,080; mg/L – 100, 150, 250, 400; ‰ - 0,200, 0,300, 0,500, 0,800).
 - In order to save, leave the screen displaying desired alarm point and wait until device emits sound signal – change has been successfully saved

Specifications

Electrochemical sensor (Fuel cell sensor)

4 digits display

LCD with backlight

Blowing time about 6 seconds

Blowing flow about 10 L/min

Blowing flow interrupt detection

Display range: 0.00~3.00‰(1.50 mg/L)

Accuracy: +/- 0.05‰ at 0.50‰

Warm-up time: 10 seconds

Testing time: 10 seconds

Repeat test: less than 60 seconds

Selective unit of measure: ‰, ‰, mg/L

Show „HI” if above the display range

LCD change color if above alcohol limit

Selective alarm point for 0.02~0.08% (default is 0.05%)

Auto power off

Test counter

Calibration interval: 1 year or 200 tests

Low battery indication

5 pcs mouthpiece included

Power supply: 1x AAA battery

Operation temp: 10°C ~ 40°C

Storage temp: 0°C ~ 50°C

Calibration and Record

- You should have your alcohol tester calibrated every 12 months or 200 tested times. Failure to do so may cause your readings to be inaccurate

Callibration Record			
Date of first use		Serial number	
#	Calibration date	Remarks	Recalibration by

Precautions

- Your blood alcohol content (BAC) may change some minutes later.
- You should have your alcohol tester calibrated every 12 months or 200 tests. Failure to do so may cause your readings to be inaccurate.
- Do not use this device to determine if you can drive.
- You cannot use the test results of device in court of law.
- You cannot use the test result of device against police's unit

Specifications are subject to change without prior notice.

English
Correct Disposal of This Product
(Waste Electrical & Electronic
Equipment)

(Applicable in the European Union and other European countries with separate collection systems) This marking shown on the product or its literature, indicates that it should not be disposed with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources. Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take this item for environmentally safe recycling. Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

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